

# DUAL LASER PORTABLE PXP-20HF & PXP-40HF



- Best choices for perfect image presentation
- Dual laser & preset focus at 26" SID
- Ultra-Light Weight & Higher Power output
- Preset APR & Remote Control Software

The **PXP-20HF PLUS** is newly designed with evolutionary engineering for equine veterinarians. Weighing only 22 lbs and a dual laser focusing on the center of x-ray film at 26" preset SID, **PXP-20HF PLUS** makes the most difficult equine radiographic imaging easier by only taking one exposure without retaking a film repeatedly.

Preset APR technique, collimator lamp and laser pointer are turned on and selected by remote hand-switch prior to making an exposure. Variable mA, 15 to 20 mA and 40 to 100 kV provide plenty of power output for penetration of shoulders, necks and C-spines.

The **PXP-40HF** is the most compact, ultra-weight and most powerful units available today. It is integrated with a single laser, but the dual laser collimator is an option. At only 27.0 lbs it is light weight and compact. **PXP-40HF** is more versatile for a small, large animal, equine or mixed animal practitioner. This unit has a range of 40-100 kV and variable 16 mA to 35 mA. It is equipped with remote control software, PROM and remote control hand-switch as standard features.

CALL OR WRITE FOR DETAILED INFORMATION  
VISIT OUR WEB SITE AT [www.UNITED-RADIOLOGY.com](http://www.UNITED-RADIOLOGY.com)



**United** Radiology Systems, Inc.

151 S. Pfingsten Rd. Unit T • Deerfield, IL 60015  
Telephone: (847)291-0999 • Fax: (847)291-1090  
E-mail: [info@UNITED-RADIOLOGY.com](mailto:info@UNITED-RADIOLOGY.com)  
Web: [www.UNITED-RADIOLOGY.com](http://www.UNITED-RADIOLOGY.com)