

PX-110CL

Battery Mobile X-Ray Unit

SPECIFICATION

Generator	Inverter control system	
Output	80 kV / 80 mA (max.) 110 kV / 60 mA (max.)	
MAs	0.3 ~100mAs [31 steps]	
X-ray tube	Focal spot	1.4mm IEC
	Anode	Stationary
	Heat unit	50 KHU
	Al. eq.	0.6mmAl.eq.Min
Buttery	Type	HV5-12
	Capacity	20 each in series
	Power	DC. 240V
	Charging time	Approx. 8 Hours
	Charge	Single phase 110V 50/60 Hz
Collimator	Manual type	
	20 sec. Timer	
Tube Arm	Gas Spring balanced	
Cassette Box	14"x17" cassette 10 pieces Option: 14"x17" cassette 15 pieces	
Shipping Guide	Packing bulky	1.06 x 0.66 x 1.55 (LxWxH)
	Net weight	150kg
	Gross weight	230kg (depends on packing materials)



UNITED RADIOLOGY SYSTEMS,INC

151 S. Pfingsten Rd. Suite T
Deerfield, IL 60015

Tel: 847-291-0999 Fax: 847-291-1090

e-mail : sales@united-radiology.com web: www.united-radiology.com